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
Original

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Instructions

D

E

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
Instructions

E

We are pleased that you decided to buy a highly developed piece of equipment from SCHICK and would like to wish you every success when working with your new Air puffer.

We wrote these instructions to enable you to get accustomed to your new piece of equipment and to provide you with the correct operating and maintenance instructions.

Index	page
1. Safety instructions.....	3
2. Range of applications.....	3
3. Summary - List of contents/Accessories - Optionals.....	4
4. Initial start up and installation.....	5, 6
5. Installation samples.....	6
6. Technical data.....	7
7. Declaration of conformity.....	8

Schick Dental presents another innovative product for the optimal workbench - the Air puffer, which is attachable directly to the suction funnel. The foot switch operated Air puffer is located directly inside the work field thus, the inconvenient taking of the standard air nozzle is not necessary any more. The limberly joints of the air puffer offer high flexibility (extensible) and the adjustment of the air flow can be done easily by the adjustment screw on the air puffer body.

In combination with the illuminated protective pane L Protect, the Air puffer is the ideal solution for every work bench to offer the technician as much comfort, efficiency and safety as possible.

1. Safety instructions

1.1 Ascertain that your mains compressed air supply coincides with the data on the technical description.

1.2 The Air puffer units are not suitable for the following applications:

- for medical applications

Caution!

- Repairs and other technical procedures must only be carried out by suitable qualified personnel, authorized by SCHICK.

- SCHICK do not guarantee the Air puffer unit should it not have been used in accordance with the instructions.

2. Range of applications

The Air puffer is made for the universal use in dental laboratories on the technicians work bench for puffing on or cooling dental materials, tools or work pieces.

Conditions of environment:

- interior 5° - 40° C.
- up to 2,000 meter over sea level

Grade of pollution: 2

WEEE-Reg.-Nr. DE 78620387

3. Summary - List of contents

(1)

(2)

(3)

(4)

Air puffer complete	Art. No.: 4400
(1) connection block	Art. No.: 4401
(2) pneumatic foot switch	Art. No.: 4412
(3) set of hoses	Art. No.: 4410
(4) connector set	Art. No.: 4413

4. Initial start up and installation

4.1 Installation of the air puffer

The Air puffer can be positioned almost everywhere at the work place due to the selfadhesive connectors stripes (4) (pic.1).

The surfaces have to be cleaned carefully! After removing the protective foil, the connector stripe has to be pressed firmly to the surfaces. When connecting the connector stripes to each other, they have to lock audible for a save connection.

pic. 1


pic. 2

4.2 Installation of the air hoses

Footswitch and Air puffer are equipped with quick release plugs for an easy and quick connecting and disconnecting of the hoses.

Scheme:


Connect air hose 1 as first to the plug 1 on the foot switch as seen on the scheme above. Afterwards connect the Air puffer to the foot switch by using hose 2 . Make sure that the hoses are plugged in completely! Now the compressed air plug NW 7,2 can be connected to the compressed air supply.

Always connect the open hose side as first to the foot switch before connecting the hose to the compressed air supply!

4.3 Adjusting the air pressure

The air pressure can be adjusted directly on the air puffer independent from the intake pressure.

For increasing or decreasing the air pressure, the rear adjusting screw (1) can be turned left- or rightwards. The adjusting screw can be locked by screw (2).

4.4 Using the footswitch

The Air puffer can be activated by the foot switch. When pushing the foot switch, air streams as long as the foot switch is pushed. It is not possible to adjust the pressure by the foot switch.

4.5 Disconnecting the hoses

The Air puffer can be separated easily from the hoses by the quick release connectors.

For disconnecting the hoses, push the ring into the connector and pull the hose as seen on the picture above.

5. Installation samples


pic. 1

pic. 2

pic. 3

pic. 4

We recommend to fix the Air puffer on the side of the suction funnel as shown on pic.1 - 4 .
Alternatively the Air puffer can be fixed directly on the work bench. Due to the joint hose the nozzle can be moved into almost every position.

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declare herewith, that the product

Air puffer 4400
Pneumatic foot switch 4412

is in conformity with the following provisions of Directive:

2006/42/EG (machinery directive)

Name and address of person in charge:	Wolfgang Schick Lehenkreuzweg 12 88433 Schemmerhofen
---------------------------------------	--

Schemmerhofen, February 2016

A handwritten signature in black ink, appearing to read 'W. Schick'.

W. Schick
General manager

Subject to technical modifications without prior notice

These instructions should be readily accessible and are best kept close to the unit itself.

We would like to take this opportunity to advise you that a proper repair service and suitable qualified personnel are required for such highly developed technical equipment. SCHICK guarantees to carry out perfect repairs using original spare parts

Schick GmbH
Lehenkreuzweg 12
D-88433 Schemmerhofen
Telefon +49 7356 9500-0
Telefax +49 7356 9500-95
E-Mail info@schick-dental.de
Internet www.schick-dental.de

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We, the SCHICK GmbH