


# C3-Master (1.Vers.)



Ersatzteilliste für Motorhandstück Art.Nr. 7000  
List of spare parts for motorhandpiece art.no. 7000



4125 Stift / pin




7102 Ø2,35mm Spannzange / chuck  
(7103 Ø3,00mm) Spannzange / chuck



7113 Spitze / tip



7118 Ø2,35mm Dichtkappe / sealing cap  
(7118/1 Ø3,00mm) Dichtkappe / sealing cap




7110 Griffhülse / grip sleeve




7144 Schraube / screw


Werkzeuge für Spannzangenwechsel  
Tools to change chuck




6223 Gegenhalteschlüssel /  
counterstay wrench



4115 Spannzangenschlüssel / chuck key



5158 Gabelschlüssel / wrench



7250 Motorkabel / motor cable



**SCHICK**  
DENTAL

Schick GmbH  
Lehenkreuzweg 12  
D-88433 Schemmerhofen  
Telefon +49 7356 9500-0  
Telefax +49 7356 9500-95  
E-Mail [info@schick-dental.de](mailto:info@schick-dental.de)  
Internet [www.schick-dental.de](http://www.schick-dental.de)